

## ImmunoPrism™ Sequencing Index Guide

ImmunoPrism Index#	Index Primer Sequence Read
1	AGTCAA
2	CCGTCC
3	GTCCGC
4	GTGGCC
5	GAGTGG
6	ACTGAT
7	CGTACG
8	ATTCCT
9	ATCACG
10	CGATGT
11	TTAGGC
12	TGACCA
13	ACAGTG
14	GCCAAT
15	CAGATC
16	ACTTGA
17	GATCAG
18	TAGCTT
19	GGCTAC
20	CTTGTA
21	AGTTCC
22	ATGTCA
23	GTGAAA
24	GTTTCG

**Note:** For low-level multiplexing for sequencing, it is recommended that the following index primers be used:  
 Pool of 2 samples: Index 3 and 6 | Pool of 3 samples: Index 1, 3 and 5 | Pool of 4 samples: Index 2, 3, 4 and 5